

### GANTRY ROTATION vs. SPECIMEN ROTATION: IMPLICATIONS FOR IN VIVO IMAGING

#### Introduction

In computed tomography (CT) and related imaging modalities, projection data for reconstruction can be acquired by either rotating the gantry around a stationary subject or rotating the subject in front of a fixed source–detector pair. From a purely geometric and mathematical standpoint, these two approaches are equivalent—assuming a rigid object and a perfectly calibrated system.

However, when applied to in vivo imaging of living animals, which are inherently non-rigid and subject to motion and deformation, the equivalence breaks down. Both published studies and current market systems indicate that biological motion, mechanical stability, and practical considerations all favor gantry rotation as the method of choice. Along with patient comfort, this is the key reason why clinical linear accelerators only employ gantry rotation.

#### System Design and Application Domain

Most in vivo small-animal CT systems employ a rotating gantry with a stationary bed. This configuration:

- Minimizes physical stress on the animal.
- Maintains stable anesthesia, monitoring, and catheterization lines.
- Reduces motion artefacts caused by physiological processes, like respiration or cardiac activity.
- Simplifies monitoring during scans, ensuring uninterrupted access to support systems (IV lines, ECG, breathing, body temperature control).
- Reduces the impact of X-ray scatter on image quality, particularly along the long axis of the animal, improving contrast and spatial resolution in projection data.

During in vivo imaging, multiple support systems are required, such as intravenous lines, respiratory assistance, and physiological monitoring (ECG, breathing, body temperature). By keeping the animal on a stationary bed, these connections remain stable and accessible, preventing torsion of cables and ensuring consistent anesthesia and monitoring throughout the scan.

#### Mechanical Stability and Calibration

Rotating specimen stages demand precise alignment of the rotational axis to minimize eccentricity and wobble, both of which may generate ring or streak artifacts. Additionally, rotating beds require motors and transmission systems that can introduce micro-vibrations or resonances within the animal.

In contrast, rotating gantries primarily require mechanical stability of the source–detector assembly, which is more easily maintained through routine calibration. Gantry rotation, particularly when physically separated from the bed, as in the SmART+ (Small Animal RadioTherapy) preclinical irradiator (Precision X-Ray, Inc., Madison, CT), causes minimal disturbance and reduces mechanical vibrations in the animal, consequently improving signal-to-noise ratio (SNR) reconstruction quality.

#### Reduction of Deformation and Motion Artifacts

Animal rotation can introduce soft-tissue deformation through centrifugal forces, imperfect fixation, or micro-vibrations of the stage. Such effects can manifest as blurring, rings, or streak artifacts in reconstructions. By keeping the subject stationary, these sources of distortion are largely avoided, leading to more reliable anatomical and functional imaging.

## Practical Considerations

For live animal studies, rotating the gantry is safer and more practical:

- Prevents torsion or tangling of anesthesia and monitoring lines.
- Reduces risk of stress or injury to the subject.

## Clinical Relevance

The choice between gantry and specimen rotation is not only important in preclinical research but also highly relevant to clinical practice. In human CT imaging, gantry rotation is universally adopted for many of the same reasons observed in animal studies:

**Translational continuity:** Using gantry-based imaging in preclinical systems mirrors the clinical standard, ensuring that animal imaging protocols, artefact characteristics, and reconstruction strategies closely align with those used in patient care. This strengthens the translational bridge between preclinical findings and clinical applications.

**Workflow efficiency:** Gantry rotation supports faster setup, easier positioning, and integration of multimodal imaging—all of which are essential in both preclinical and clinical environments.


**Minimized motion artefacts:** Keeping the patient stationary reduces the possibility of misalignment during image capture, which is critical for high-quality diagnostic imaging.

**Patient comfort and safety:** Just as small animals require stable anesthesia lines, human patients require stable IV access, ventilators, ECG leads, and monitoring devices during imaging. A rotating bed could compromise these connections and pose risks.

By adopting gantry rotation in small-animal imaging, researchers not only optimize image quality and animal welfare but also ensure methodological continuity with clinical imaging practices, enhancing the translational impact of their work.

## Summary

While gantry and specimen rotation may be mathematically equivalent in theory, they are not interchangeable in practice for *in vivo* imaging. Gantry rotation provides superior stability, reduces artifacts, and ensures animal safety by minimizing stress and preserving physiological monitoring. For these reasons, rotating gantry designs are the established standard in both commercial small-animal imaging systems and published preclinical studies.



**Figure 1.** ROTATING GANTRY SmART+ (Small Animal Radiotherapy System), Precision X-Ray, Inc. (Madison, CT, USA). On the left, the gantry (orange) is shown with the X-ray tube positioned above and the detector below the specimen at the initial 0° angle. On the right, the gantry rotates around the specimen, which remains immobile at the center.

## References

- Dillenseger, J. P., Goetz, C., Sayeh, A., Healy, C., Duluc, I., Freund, J. N., ... & Choquet, P. (2017). Estimation of subject coregistration errors during multimodal preclinical imaging using separate instruments: origins and avoidance of artifacts. *Journal of Medical Imaging*, 4(3), 035503-035503.
- Hsieh, J. (2024). Spatial and temporal motion characterization for x-ray CT. *Medical Physics*, 51(7), 4607-4621.
- Kyme, A. Z., & Fulton, R. R. (2021). Motion estimation and correction in SPECT, PET and CT. *Physics in Medicine & Biology*, 66(18), 18TR02.
- Nardi, C., Molteni, R., Lorini, C., Taliani, G. G., Matteuzzi, B., Mazzoni, E., & Colagrande, S. (2016). Motion artefacts in cone beam CT: an in vitro study about the effects on the images. *The British journal of radiology*, 89(1058), 20150687.